
pulverisette 15



Cutting Mill

- Suitable for soft to medium-hard dry materials
- Rapid, reproducible grinding
- Easy to operate and safe to use

grinding sample cutting mill preparation for your lab



Cutting Mill

Field of application

Grinding of soft to medium-hard dry materials and fibrous materials or materials containing cellulose including rubber, leather, paper, cardboard, textiles, plastics, wood, metal-free waste, coal, feedstuffs, maize, grain, confectionery, malt, noodles, herbs and spices, dried meat, bone, horn, pills, tablets, foliage, fibres, peat, roots, tobacco, cork.

Straw, stalks foils and similar long, thin samples can be fed full-length into the instrument.

Method of operation

In the grinding chamber of the pulverisette 15 cutting mill, four knives rotate against three stationary knives in the cutting action to reduce the sample. The material remains in the grinding chamber until the desired final particle size is achieved and the continuous flow of air carries it through the mesh of the interchangeable sieve insert into the collecting vessel.

Advantages

- Easy cleaning of grinding chamber and grinding parts
- Convenient materials feed
- Rapid, reproducible grinding
- Defined final particle size with interchangeable sieve inserts
- Ease of maintenance
- Certified safety (CE mark)
- 2 year guarantee




quality fine grinding control Cutting Mill





pulverisette 15 with collecting vessel for large quantities



cutting rotor and sieve inserts



grinding chamber

Design Characteristics

- Easy access to grinding chamber using only two knurled screws
- Housing with cover opens out completely
- Cover can be removed from hinge
- 2-litre stainless steel funnel with plastic pusher
- Grinding chamber efficiently sealed
- Long lasting hardened steel knives
- Cutting edges can be reused
- Knives can be resharpened
- Cutting gap externally adjusted
- High cutting rotor speed (2800 rpm)
- Bearing support on both sides of cutting rotor
- Dust-sealed bearings
- Recyclable, torsion-resistant aluminium housing
- Designed for operation either fixed to bench-top or for use on independent mounting frame.

Accessories

■ Sieve inserts

Sieve insert made of stainless steel with trapezoidal or square perforations can be selected from the range to obtain final particle size of your choice.

For free-standing installation of the pulverisette 15

Filter hose for large quantities

A 30 litre vessel connected by a filter hose provides

www. milling A 30 litre vessel connected bulk collection facility.

pulverisette 15

fritsch.de

pulverisette 15

Technical data

| min. sample quantity | 30 - 50 ml |
|--|--|
| max. throughput (depending on the material and sieve size) | 50 l/h |
| final fineness depending on the selected sieve insert | 0.25 - 6 mm |
| feeding | batchwise/continuous |
| rotor speed | 2,800 - 3,400 rpm (depending on voltage and frequency) |
| electrical details | 400 V/3~, 50 Hz, 1900 watt 230-240 V/1~, 50 Hz, 2100 watt 100-120 V/1~, 60 Hz, 1800 watt |
| motor-shaft-power according to VDE 0530, EN 60034 | 1.1 kW for 100-120 V/1~ motor 1.5 kW for all other motors |
| weight | net: 42 kg, gross: 72 kg |
| dimensions w x d x h | table mounting or on stand: 42 x 48 x 69 cm |
| packing details | wooden case 87 x 58 x 78 cm |

Ordering data

| Order no. | Description |
|-------------|--|
| | Cutting Mill pulverisette 15 without sieve insert, incl. 3.5 I collecting vessel * included in the basic price of the instrument |
| 15.4030.00 | for 400 V/3~, 50 Hz, 1900 watt Attention: The pulverisette 15 with voltage "400 V/3~" can <u>only</u> be operated on a three phase supply network! |
| 15.4020.00 | for 230-240 V/1~, 50 Hz, 2100 watt |
| 15.4010.00 | for 100-120 V/1~, 50-60 Hz, 1800 watt |
| | Other voltages on request! |
| | Sieve inserts made of stainless steel |
| 45.1200.10 | 0.25 mm trapezoidal perforation |
| 45.1210.10 | 0.5 mm trapezoidal perforation |
| 45.1220.10 | 0.75 mm trapezoidal perforation |
| 45.1230.10 | 1 mm trapezoidal perforation |
| 45.1240.10 | 1.5 mm trapezoidal perforation |
| 45.1250.10 | 2 mm trapezoidal perforation |
| 45.2200.10 | 4 mm square perforation |
| 45.2210.10 | 6 mm square perforation |
| | Collecting vessels |
| 15.4400.00* | 3.5 litres |
| 15.4560.00 | 60 litres with filter hose |
| | Accessories/Spare parts |
| 83.4300.00 | stand |
| 45.4000.00* | set rotating knives made of tempered steel (set = 4 pieces) |
| 45.4100.00* | set stationary knives made of tempered steel (set = 3 pieces) |



Fritsch GmbH Manufacturers of Laboratory Instruments Industriestrasse $8\cdot55743$ Idar-Oberstein \cdot Germany

Phone: +49 67 84 70 0 \cdot Fax: +49 67 84 70 11 E-Mail: info@fritsch.de \cdot Internet: www.fritsch.de